

# Steps for instrument Cleaning and Disinfecting

1

Soak



Soak used instruments

2

Clean



Clean using detergent, gloves and a brush

3


Rinse



Rinse with clean water

4

Dry



Air dry or dry with lint free cloth

5

Immerse



Immerse in approved disinfectant. Follow manufacturers instructions

6


Rinse



Rinse with clean water

7

Dry



Air dry or dry with lint free cloth

8

Store



Store in a dry, clean space



**simcoe  
muskoka**  
DISTRICT HEALTH UNIT

Tel: 705-721-7520  
Toll free: 1-877-721-7520  
[www.simcoemuskokahealth.org](http://www.simcoemuskokahealth.org)  
*Your Health Connection*